

Signal Enhancement in Magnetic Particle Imaging by a passive Dual Coil Resonator (pDCR)

Bachelor Thesis Defence Talk

Jonas Philipps

16.02.2021



1 Basics and Motivation

2 Principles of the pDCR

3 Impact of the pDCR on SNR and Spatial Resolution


4 Conclusion & Outlook

- Novel tracer-based imaging technique
- Determination of the spatial distribution of a tracer by stimulation with magnetic fields
- Tracer: Material to be imaged
 - Superparamagnetic iron oxide particles (SPIONs)
- Advantages:
 - High temporal resolution → Real-time imaging
 - High sensitivity
 - Good spatial resolution
 - Quantitative
 - No use of harmful ionizing radiation




SPIONs

- Superparamagnetism → ferromagnetism + single domain
⇒ no particle-particle interactions
- Nonlinear magnetization behavior → Langevin function




Signal Generation

- Excitation with a sinusoidal magnetic field (drive field)



- Detection via voltage induced into receive coils



Spatial Encoding

- So far: Magnetization response of all particles
→ Mean concentration
- Wanted: Particle distribution
→ Local concentration
- Solution: Superposition of a gradient field (selection field) featuring a field free point (FFP)






Image Reconstruction

- $\mathbf{Sc} = \hat{\mathbf{u}}$
 - \mathbf{c} : particle concentration
 - $\hat{\mathbf{u}}$: induced voltage
 - \mathbf{S} : System matrix
 - Describes the spectra of a delta sample being placed at every position in the volume to be imaged
- In true measurements: $\tilde{\mathbf{u}}$ distorted by noise
 - $\|\mathbf{Sc} - \tilde{\mathbf{u}}\|_2^2 \xrightarrow{c} \min$
 - Solve with an iterative solver (Kaczmarz method)
- \mathbf{S} also noisy
 - Denoising: Fourier or cosine transformation will compress the true signal but not the noise
→ Noise can be separated by thresholding

■ Scanner:

- 1) selection field coils (blue): generate the FFP
- 2) drive field coils (green): generate the excitation field
- 3) receive coils (red): pick up the particles' magnetization response


■ Signal chain:




[4]

Motivation – Signal Enhancement


- True measurements: Noisy signal
→ Higher harmonics disappear in a noise floor




- More harmonics detectable
→ Better image reconstruction possible
- Spatial information is encoded in higher harmonics
→ Exemplarily shown by the Fourier series of the square wave function




(a) Original square wave function



(b) Fourier Series with $n=5$ elements




(c) Fourier Series with $n=15$ elements




(d) Fourier Series with $n=100$ elements

Basics – The passive Dual Coil Resonator (pDCR)

- Goal: Enhancement of SNR at a high frequency band
→ improve spatial resolution
- Implementation: Passive coaxial dual coil resonator
 - Inner coil (50 mm) to receive SPION signal with high sensitivity
 - Outer coil (110 mm) coupled into native receive coil → combined inductance $210\ \mu\text{H}$
 - Coils connected via a capacitor ($0.314\ \text{pF}$)
→ resonance at 608 kHz



(b)



(c)




- System matrix simulation:

$$S_{ij} = i\omega_i \mu_0 m \mathbf{R} \frac{V}{N} \mathcal{F}[\mathcal{L}(\mathbf{r}(t) - \mathbf{x}_j)](\omega_i) \cdot \hat{\Gamma}(\omega_i)$$


- Signal simulation:

- (a) Magnetization signal in frequency domain
- (b) Magnetization signal in time domain
- (c) Voltage signal in time domain
- (d) Voltage signal in frequency domain



■ Transfer function:


- Transfer function of receive chain incl. pDCR
 - Increased sensitivity at high frequency band
 - Peak at 694 kHz
 - inside scanner: coupling with scanner elements
- resonance frequency shifts





■ Signal-to-noise ratio:

$$SNR = \frac{S(f) - BG(f)}{\sigma(BG(f))}$$

- SNR of system matrix used for threshold determination in image reconstruction





- Threshold variation
- Determination of number of harmonics with amplitude above threshold
 - (a) pDCR gives advantage for all thresholds
 - (b) Ratio $n_h(\text{pDCR})/n_h(\text{native})$ peaks at a threshold of 0.8





Impact of the pDCR – Threshold and Image Reconstruction

Threshold: 0.3





- Low threshold: No difference, both reconstructions perform well

Threshold: 2.5



- Medium threshold: pDCR reconstruction performs well, while native gets worse


Threshold: 4.0





- High threshold: pDCR reconstruction gets worse, but still better than native

Impact of the pDCR – Shift of Resonance Frequency



- Goal: Increase the total number of harmonics above threshold
 - Shifting the resonance frequency of the pDCR to higher frequencies
 - Results:
 - Higher resonance frequency
 - Higher maximum ratio
 - Lower threshold at which the maximum is reached
- ⇒ Shifting the resonance frequency can optimize the performance of the pDCR
- ⇒ It depends on the threshold, which resonance frequency performs best




- 3D drive field
- Only signal induced in x -channel used
- Signal simulation:
 - (a) Magnetization signal in frequency domain
 - (b) Magnetization signal in time domain
 - (c) Voltage signal in time domain
 - (d) Voltage signal in frequency domain
- Transfer function of receive chain for signal tuning



- Threshold variation
- Determination of number of harmonics with amplitude above threshold
 - (a) pDCR gives advantage for all thresholds
 - (b) Ratio $n_h(\text{pDCR})/n_h(\text{native})$ peaks at 0.9




- Image reconstruction of a delta sample with a threshold of 0.9





- Reconstructed intensity increased
- Signal flattens faster
⇒ Spatial resolution has improved

- Goal: Increase the total number of harmonics above threshold
- Shifting the resonance frequency of the pDCR to higher frequencies
- Results:
 - ⇒ Shifting the resonance frequency can optimize the performance of the pDCR
 - ⇒ It depends on the threshold, which resonance frequency performs best



- **Scanner:** Bruker Preclinical MPI system
 - Drive field amplitude: 14 mT
 - Drive field frequencies: 24.5 kHz, 26 kHz, 25.3 kHz
 - Max. magnetic field gradient: 2.5 Tm^{-1}
- **Tracer:**
 - C2 (in-house, hydrod. diameter: 49.4 nm, iron concentration: 8.5 mg/ml)
 - Perimag (commercial, hydrod. diameter: 130 nm, iron concentration: 8.5 mg/ml)
- **Measurements:**
 - SNR measurements
 - (a) C2
 - (b) Perimag
 - Resolution measurements





Impact of the pDCR – Resolution Phantom

MPI Resolution phantom

Material: Acrylic glass (unpolished)

All dimensions in mm




Impact – Reconstructed Images (C2, ParaVision)

- Signal normalized
→ No comparison of reconstructed intensities possible
- No system matrix denoising → 500 averages for SM


⇒ Spatial resolution improved by the pDCR

- Points in clusters can be better resolved
- Individual points have an increased intensity compared to a cluster (mainly observable in the prominent plane)


(a) Native – prominent plane




(b) pDCR – prominent plane



(c) Native - overlay



(d) pDCR - overlay




- Signal normalized
→ No comparison of reconstructed intensities possible
- System matrix denoising → number of averages reduced to 300


⇒ Spatial resolution improved by the pDCR

- Points in clusters can be better resolved
- Individual points have an increased intensity compared to a cluster (mainly observable in the prominent plane)


(a) Native – prominent plane




(b) pDCR – prominent plane



(c) Native - overlay



(d) pDCR - overlay




- Signal not normalized
→ Comparison of reconstructed intensities possible
- System matrix denoising


⇒ Spatial resolution improved by the pDCR

- Points in clusters can be better resolved
- Individual points have an increased intensity (mainly observable in the prominent plane)


(a) Native – prominent plane




(b) pDCR – prominent plane



(c) Native - overlay



(d) pDCR - overlay




- Signal not normalized
→ Comparison of reconstructed intensities possible
- System matrix denoising


⇒ Spatial resolution improved by the pDCR

- Points in clusters can be better resolved
- Individual points have an increased intensity (also observable in the overlay)


(a) Native – prominent plane




(b) pDCR – prominent plane



(c) Native - overlay



(d) pDCR - overlay



■ Conclusion:

- Simulation (1D and 3D drive field):
 - ✓ pDCR enhances the SNR of higher harmonics
 - ✓ FWHM of PSF (spatial resolution) is increased by pDCR
- Shifting the resonance frequency to higher harmonics:
 - ✓ Number of harmonics above threshold increases
 - ✓ Ratio $n_h(\text{pDCR})/n_h(\text{native})$ depends on threshold
- Measurements (at Bruker scanner):
 - ✓ SNR enhancement of higher harmonics validated (Perimag and C2)
 - ✓ Resolution phantom: spatial resolution and absolute intensities increased with the pDCR (Perimag and C2)

■ Outlook:

- Optimizing noise modelling (especially in 3D simulation)
- Using all three channels in 3D simulation
- Different pDCR designs:
 - Tunable pDCR
 - 3D pDCR
 - Mini pDCR

Thank you for your attention!

Any questions?


RWTHAACHEN
UNIVERSITY

ExMI Experimental
Molecular
Imaging

PMI PHYSICS OF
MOLECULAR
IMAGING SYSTEMS

- [1]** Tobias Knopp and Thorsten M. Buzug. Magnetic Particle Imaging: An Introduction to Imaging Principles and Scanner Instrumentation. Heidelberg: Springer, 2012
- [2]** L. W. E. Starmans & H. Gruell, PLOS, 2013
- [3]** Kerstin Lüdtke-Buzug, University of Lübeck
- [4]** T Knopp, N Gdaniec, and M Möddel. “Magnetic Particle Imaging: From Proof of Principle to Preclinical Applications”
- [5]** S D Reinartz et al. “Feasibility of a Spatial Resolution Enhancement by a Passive Dual Coil Resonator (pDCR) Insert for Large Bore MPI Systems”

Selection Field and Nonlinear Magnetization Behavior



V-Phantom Carrier

